



Dornier **COMPACT SIGMA**

A truly modular lithotripter



Why Compact Sigma



When combined with other devices, such as C-arms, ultrasound scanners, and patient tables, this device transforms swiftly into a multi-functional urological solution.




Mobile and compact lithotripter for varied treatment needs



Excellent system efficiency for optimized treatment outcomes



Various imaging options for your ideal solution



Mobile and compact lithotripter for varied treatment needs



Remarkable capabilities

- Enhance lithotripsy practice through customized configuration to fit your unique needs
- QuickLinx – Dornier’s integrated optical position control, for ease of C-arm alignment and faster procedure preparation

Ergonomic design

- Easily transportable with large wheels and foldable therapy head
- Small footprint; requires minimal space in any operating room or storage facility

Plug and play

- No installation required; plug in and start operating instantly

Wide range of applications

- Specially designed for ESWL and ESWT*

*Refer to slide 7 for a list of indications

Excellent system efficiency for optimized treatment outcomes



Dornier's EMSE technology


- Most clinically proven shock wave technology, Dornier's exclusive¹
- Offers consistent quality performance with reliable energy output

Flexible therapy head positioning

- Facilitates large window of access to various stone locations with 36 therapy head positions
- Ensures patient comfort with over- and under-table positions as patients, particularly the elderly and obese, can remain supine throughout the procedures



Various imaging options for your ideal solution



Comprehensive suite of solutions

- Customize your imaging preferences with Dornier's trusted and qualified partners
- Achieve outstanding image quality with high-end Flat Panel Detector C-arms

Ultrasound-guided ESWL

- Eliminates radiation exposure for patients and users while allowing users to detect any type of stones
- Real-time monitoring offers immediate readjustment of patient positions when necessary and maintains stone alignment over the target zone
- Isocentric ultrasound imaging provides a high degree of flexibility and imaging quality²

Various imaging options for your ideal solution



Dual imaging option: The best of both worlds

- Comprehensive imaging capabilities, with the option of localizing stones using either X-ray or ultrasound imaging, or both simultaneously
- X-ray stone localization allows fast, initial positioning of the patient, while ultrasound provides real-time information on patient respiration and stone disintegration



Offering a broad range of applications

The Dornier Compact Sigma can be used for the following indications³:

Extracorporeal Shock Wave Lithotripsy (ESWL)

- Renal pelvis and renal calyx calculi
- Upper and lower ureteral calculi
- Prevesical ureteral calculi
- Biliary calculi, excluding gall bladder and cystic duct stones
- Bile duct calculi

Extracorporeal Shock Wave Therapy (ESWT)

- Pancreatic calculi
- Impaired fracture healing and pseudarthrosis



Did you know?

The Compact Sigma can be operated remotely. This decreases the exposure to ionizing radiation for users and other healthcare personnel.



Dornier *MedTech*

References:

¹Data on file at Dornier MedTech

²Macchione, N., Elia, A., Gofrit, O., Pode, D., & Duvdevani, M. (2013). SWL with continuous targeting by ultrasound; are there benefits? European Urology Supplements, 3(12), 51. [https://doi.org/10.1016/S1569-9056\(13\)61727-3](https://doi.org/10.1016/S1569-9056(13)61727-3)

³Refer to the Dornier-Compact Sigma Operating Manual for more information. Indications may differ depending on region.

The information contained in this material has been obtained directly from the Dornier Compact Sigma® brochure DMT746-SIGMA-BROC-REVA-0323-GL-EN.

©2023 Dornier MedTech. All rights reserved. The contents herein and product shown (including its design, specifications and functions) are subject to change without prior notice. Dornier Compact Sigma® is a registered trademark of Dornier MedTech and its use throughout this document is protected. The information contained in this material is for information purposes only and provided "as is" without any warranties of any kind, whether expressed or implied, including but not limited to, implied warranties of satisfactory quality, fitness for a particular purpose and/or correctness. The actual product may vary from the images shown. The product shown is only offered in countries where it may be legally offered. Please check the availability with your local sales representative or customer service. DMT746-SIGMA-PPT-REVA-0323-GL-EN